

Vertebral Fracture Assessment

John Schousboe, MD, PhD
Park Nicollet Clinic & Specialty Center

The current generation of Hologic DXA scanners support low dose single and dual energy x-ray imaging of the AP and lateral spine to allow detection of vertebral fractures, including those that were not recognized clinically at the time of their occurrence. Originally implemented on Hologic fan beam densitometers as the HD Instant Vertebral Fracture™ assessment (IVA HD), in 2006 the International Society of Clinical Densitometry (ISCD) recommended that the term “Vertebral Fracture Assessment” (VFA) be used to refer to this imaging modality.¹

Vertebral fractures have long been recognized as a cardinal sign of systemic skeletal fragility and predictor of subsequent vertebral and non-vertebral osteoporotic fractures, even after accounting for age and bone mineral density (BMD).²⁻⁵ Recent US population-based data using the US National Health Nutrition Examination and Survey (NHANES) showed the prevalence of vertebral fractures rises sharply with increasing age from 2.1% among those age 40 to 49 to 18.0% among those age 80 and older, with the prevalence for men and women within each age group being nearly the same.⁶ However, other studies suggest that the age-adjusted prevalence of vertebral fractures may be higher in older women compared to older men.^{5,7-9}

Both the National Osteoporosis Foundation (NOF)¹⁰ and the ISCD¹¹ have established guidelines for VFA to identify patients with a prevalent vertebral fracture. Over the last 20 years, the scientific evidence supporting use of VFA at the point of care of bone densitometry to identify those with this important fracture risk factor has been well established. First, numerous studies have shown that patients with prevalent vertebral fractures (especially those that are moderate or severe) can be identified accurately on VFA images compared to standard radiography.¹²⁻¹⁴ Second, NOF and ISCD guidelines on VFA select a population for imaging with a reasonable probability of a prevalent vertebral fracture. Studies have shown by following NOF or ISCD guidelines, clinicians can capture at least half of those in the overall population who have a prevalent vertebral fracture.^{8,10,15,16} Third, cohort studies have shown that those who have

a vertebral fracture specifically identified on VFA are at higher risk of: subsequent vertebral,^{17,18} hip,^{17,19} major osteoporotic,¹⁷ non-vertebral,¹⁷ all clinical fractures,¹⁹ and fracture hospitalizations,¹⁸ even after accounting for BMD, age, and other risk factors.

VFA has several significant advantages compared to standard radiographs. First, VFA exposes patients to only 1/50th of the radiation exposure of standard thoracic and lumbar spine radiographs. Second, the fan beam technology allows visualization of vertebrae with less parallax than standard cone-beam radiography. A disadvantage of VFA is that it can be more difficult to visualize vertebral fractures superior to T7 because of the lower spatial resolution of VFA compared to standard radiography. However, this disadvantage is mitigated by the fact that osteoporotic vertebral fractures occur primarily in the skeletal region of T10 to L2, and that vertebral fractures above T7 are uncommon.^{20,21}

The third and perhaps most important advantage of VFA is that it can be seamlessly and cost-effectively²² incorporated into standard bone densitometry screening routines with minimal additional time and expense. By following NOF and ISCD guidelines, bone density technologists have demonstrated the ability to accurately identify those who meet the criteria to have a VFA as part of their bone density test,^{23,24} sparing the physician ordering the bone density test the need to determine *a priori* if a VFA should be performed. Non-radiologist physicians can easily be trained in a half day course (available through both the International Osteoporosis Foundation (IOF) and the ISCD) to accurately read VFA images in a manner similar to radiologists. When criteria very similar to the NOF or ISCD guidelines are applied, studies across four different health care delivery organizations in the US and Canada showed that about 20% of those who have a bone density test will meet the criteria for a VFA and will have the procedure done within the context of actual clinical practice.^{17,24} Among these, the prevalence of vertebral fracture has been shown to be in the range of 13% to 21% at these institutions.

In addition, identifying a prevalent vertebral fracture on VFA in clinical practice has been shown to increase appropriate prescription of fracture prevention medication by physicians and use of such medication by patients.²³ The influence of VFA on medication use behavior is particularly impactful among those who have BMD in the osteopenic range, or who have a FRAX[®] estimated risk of fracture (before vertebral fracture status is known) less than 20%.²³ Even among patients who have osteoporosis by BMD criteria, identification of prevalent vertebral fractures may increase patients' perceived need for fracture prevention medication.²³ Hence, VFA can play a significant role in closing the osteoporosis treatment gap by identifying those at high risk of fracture and motivating patient compliance with medication.²³

Use of VFA to identify abdominal aortic calcification

Abdominal aortic calcification (AAC) identified on lateral spine x-rays is an indicator of systemic atherosclerosis, and has been shown in the Framingham and Rotterdam studies to be associated with subsequent cardiovascular disease events such as: coronary heart disease,²⁵ myocardial infarction,²⁶ stroke,^{27,28} congestive heart failure,²⁹ symptomatic peripheral vascular disease,²⁸ and cardiovascular disease mortality,²⁵ even after accounting for clinical cardiovascular disease risk factors and/or Framingham Heart Score. AAC can be accurately scored on the HD Instant Vertebral Fracture™ assessment from Hologic,³⁰⁻³⁴ and several studies have shown that AAC scored on these images also predict: incident myocardial infarction and/or stroke,^{35,36} incident cardiovascular disease events,^{37,38} hospitalizations for cardiovascular disease,³⁷ and all-cause mortality.³⁷

AAC has also been shown to be a risk factor in most³⁹⁻⁴³ (but not all)^{44,45} studies for hip and non-spine fractures for the subsequent five years after its identification, even after accounting for other fracture risk factors such as age, BMD, and previous fractures. Further studies are planned to better define the association of AAC on VFA images and fracture outcomes.⁴⁶

While routine use of VFA imaging for cardiovascular disease risk assessment cannot be recommended at this time, there is significant potential to identify individuals with asymptomatic atherosclerosis on the HD Instant Vertebral Fracture™ assessment whose high risk for incident cardiovascular disease might otherwise not be recognized.⁴⁷ Further studies are needed on the impact of a high level of AAC identified on the HD Instant Vertebral Fracture™ assessment. Such studies should assess how AAC may impact cardiovascular disease risk management by health care providers, their patients, and the subsequent impact on cardiovascular disease outcomes. However, today it is possible to opportunistically identify those with a high level of AAC on the HD Instant Vertebral Fracture™ assessment, thus indicating a higher risk of incident cardiovascular disease and hip fracture.

Summary

VFA imaging at the time of a bone density test is a convenient and efficient way of identifying patients with one or more prevalent vertebral fractures. This aids in identifying those at high risk of subsequent fracture and for whom fracture prevention therapies are indicated. Interpreting these images for vertebral fracture is easy to learn through the Vertebral Fracture Recognition courses of the ISCD and IOF. The HD Instant Vertebral Fracture™ assessment can also be used to identify those with a high level of AAC who are at high risk for cardiovascular disease events.

References

1. Vokes T, Bachman D, Baim S, et al. Vertebral fracture assessment: the 2005 ISCD Official Positions. *Journal of clinical densitometry : the official journal of the International Society for Clinical Densitometry*. 2006;9(1):37-46.
2. Siris ES, Genant HK, Laster AJ, Chen P, Misurski DA, Kregge JH. Enhanced prediction of fracture risk combining vertebral fracture status and BMD. *Osteoporos Int*. 2007;18(6):761-770.
3. van der Klift M, de Laet CE, McCloskey EV, et al. Risk factors for incident vertebral fractures in men and women: the Rotterdam Study. *J Bone Miner Res*. 2004;19(7):1172-1180.
4. Melton LJ, 3rd, Atkinson EJ, Cooper C, O'Fallon WM, Riggs BL. Vertebral fractures predict subsequent fractures. *Osteoporos Int*. 1999;10(3):214-221.
5. Black DM, Arden NK, Palermo L, Pearson J, Cummings SR. Prevalent vertebral deformities predict hip fractures and new vertebral deformities but not wrist fractures. Study of Osteoporotic Fractures Research Group. *J Bone Miner Res*. 1999;14(5):821-828.
6. Cosman F, Kregge JH, Looker AC, et al. Spine fracture prevalence in a nationally representative sample of US women and men aged ≥ 40 years: results from the National Health and Nutrition Examination Survey (NHANES) 2013-2014. *Osteoporos Int*. 2017;28(6):1857-1866.
7. Kwok AW, Gong JS, Wang YX, et al. Prevalence and risk factors of radiographic vertebral fractures in elderly Chinese men and women: results of Mr. OS (Hong Kong) and Ms. OS (Hong Kong) studies. *Osteoporos Int*. 2013;24(3):877-885.
8. Schousboe JT, Rosen HR, Vokes TJ, et al. Prediction models of prevalent radiographic vertebral fractures among older men. *J Clin Densitom*. 2014;17(4):449-457.
9. Johansson H, Siggeirsdottir K, Harvey NC, et al. Imminent risk of fracture after fracture. *Osteoporos Int*. 2017;28(3):775-780.
10. Cosman F, de Beur SJ, LeBoff MS, et al. Clinician's Guide to Prevention and Treatment of Osteoporosis. *Osteoporos Int*. 2014;25(10):2359-2381.
11. Rosen HN, Vokes TJ, Malabanan AO, et al. The Official Positions of the International Society for Clinical Densitometry: vertebral fracture assessment. *Journal of clinical densitometry : the official journal of the International Society for Clinical Densitometry*. 2013;16(4):482-488.
12. Malgo F, Hamdy NAT, Ticheler C, et al. Value and potential limitations of vertebral fracture assessment (VFA) compared to conventional spine radiography: experience from a fracture liaison service (FLS) and a meta-analysis. *Osteoporos Int*. 2017;28(10):2955-2965.
13. Fuerst T, Wu C, Genant HK, et al. Evaluation of vertebral fracture assessment by dual X-ray absorptiometry in a multicenter setting. *Osteoporos Int*. 2009;20(7):1199-1205.
14. Schousboe JT, Debold CR. Reliability and accuracy of vertebral fracture assessment with densitometry compared to radiography in clinical practice. *Osteoporosis international : a journal established as result of cooperation between the European Foundation for Osteoporosis and the National Osteoporosis Foundation of the USA*. 2006;17(2):281-289.
15. Schousboe JT, Rosen HR, Vokes TJ, et al. Prediction Models of Prevalent Radiographic Vertebral Fractures Among Older Women. *J Clin Densitom*. 2014;17(3):378-385.
16. Chou SH, Vokes TJ, Ma SL, Costello M, Rosen HR, Schousboe JT. Simplified Criteria for Selecting Patients for Vertebral Fracture Assessment. *J Clin Densitom*. 2014;17(3):386-391.
17. Schousboe JT, Lix LM, Morin SN, et al. Prevalent vertebral fracture on bone density lateral spine (VFA) images in routine clinical practice predict incident fractures. *Bone*. 2019;121:72-79.
18. Prince RL, Lewis JR, Lim WH, et al. Adding Lateral Spine Imaging for Vertebral Fractures to Densitometric Screening: Improving Ascertainment of Patients at High Risk of Incident Osteoporotic Fractures. *Journal of bone and mineral research : the official journal of the American Society for Bone and Mineral Research*. 2019;34(2):282-289.
19. McCloskey EV, Vasireddy S, Threlkeld J, et al. Vertebral fracture assessment (VFA) with a densitometer predicts future fractures in elderly women unselected for osteoporosis. *J Bone Miner Res*. 2008;23(10):1561-1568.
20. Lentle BC, Berger C, Probyn L, et al. Comparative Analysis of the Radiology of Osteoporotic Vertebral Fractures in Women and Men: Cross-Sectional and Longitudinal Observations from the Canadian Multicentre Osteoporosis Study (CaMos). *J Bone Miner Res*. 2018;33(4):569-579.
21. Ferrar L, Roux C, Felsenberg D, Gluer CC, Eastell R. Association between incident and baseline vertebral fractures in European women: vertebral fracture assessment in the Osteoporosis and Ultrasound Study (OPUS). *Osteoporos Int*. 2012;23(1):59-65.
22. Schousboe JT, Ensrud KE, Nyman JA, Kane RL, Melton LJ, 3rd. Cost-effectiveness of vertebral fracture assessment to detect prevalent vertebral deformity and select postmenopausal women with a femoral neck T-score > -2.5 for alendronate therapy: a modeling study. *Journal of clinical densitometry : the official journal of the International Society for Clinical Densitometry*. 2006;9(2):133-143.
23. Schousboe JT, Lix LM, Morin SN, et al. Vertebral Fracture Assessment Increases Use of Pharmacologic Therapy for Fracture Prevention in Clinical Practice. *Journal of bone and mineral research : the official journal of the American Society for Bone and Mineral Research*. 2019;34(12):2205-2212.
24. Schousboe J, McKiernan F, Fuehrer J, Binkley N. Use of a performance algorithm improves utilization of vertebral fracture assessment in clinical practice. *Osteoporosis international : a journal established as result of cooperation between the European Foundation for Osteoporosis and the National Osteoporosis Foundation of the USA*. 2014;25(3):965-972.
25. Wilson PW, Kauppila LI, O'Donnell CJ, et al. Abdominal aortic calcific deposits are an important predictor of vascular morbidity and mortality. *Circulation*. 2001;103(11):1529-1534.
26. van der Meer IM, Bots ML, Hofman A, del Sol AI, van der Kuip DA, Witteman JC. Predictive value of noninvasive measures of atherosclerosis for incident myocardial infarction: the Rotterdam Study. *Circulation*. 2004;109(9):1089-1094.
27. Hollander M, Hak AE, Koudstaal PJ, et al. Comparison between measures of atherosclerosis and risk of stroke: the Rotterdam Study. *Stroke*. 2003;34(10):2367-2372.
28. Levitzky YS, Cupples LA, Murabito JM, et al. Prediction of intermittent claudication, ischemic stroke, and other cardiovascular disease by detection of abdominal aortic calcific deposits by plain lumbar radiographs. *Am J Cardiol*. 2008;101(3):326-331.
29. Walsh CR, Cupples LA, Levy D, et al. Abdominal aortic calcific deposits are associated with increased risk for congestive heart failure: the Framingham Heart Study. *Am Heart J*. 2002;144(4):733-739.

30. Setiawati R, Di Chio F, Rahardjo P, Nasuto M, Dimpudus FJ, Guglielmi G. Quantitative Assessment of Abdominal Aortic Calcifications Using Lateral Lumbar Radiograph, Dual-Energy X-ray Absorptiometry, and Quantitative Computed Tomography of the Spine. *J Clin Densitom.* 2016;19(2):242-249.
31. Schousboe JT, Wilson KE, Hangartner TN. Detection of aortic calcification during vertebral fracture assessment (VFA) compared to digital radiography. *PLoS ONE.* 2007;2(1):e715.
32. Schousboe JT, Wilson KE, Kiel DP. Detection of abdominal aortic calcification with lateral spine imaging using DXA. *Journal of clinical densitometry : the official journal of the International Society for Clinical Densitometry.* 2006;9(3):302-308.
33. Toussaint ND, Lau KK, Strauss BJ, Polkinghorne KR, Kerr PG. Determination and validation of aortic calcification measurement from lateral bone densitometry in dialysis patients. *Clin J Am Soc Nephrol.* 2009;4(1):119-127.
34. Cecelja M, Frost ML, Spector TD, Chowienczyk P. Abdominal aortic calcification detection using dual-energy X-ray absorptiometry: validation study in healthy women compared to computed tomography. *Calcif Tissue Int.* 2013;92(6):495-500.
35. Bolland MJ, Wang TK, van Pelt NC, et al. Abdominal aortic calcification on vertebral morphometry images predicts incident myocardial infarction. *J Bone Miner Res.* 2010;25(3):505-512.
36. Schousboe JT, Taylor BC, Kiel DP, Ensrud KE, Wilson KE, McCloskey EV. Abdominal aortic calcification detected on lateral spine images from a bone densitometer predicts incident myocardial infarction or stroke in older women. *Journal of bone and mineral research : the official journal of the American Society for Bone and Mineral Research.* 2008;23(3):409-416.
37. Lewis JR, Schousboe JT, Lim WH, et al. Long-Term Atherosclerotic Vascular Disease Risk and Prognosis in Elderly Women With Abdominal Aortic Calcification on Lateral Spine Images Captured During Bone Density Testing: A Prospective Study. *J Bone Miner Res.* 2018.
38. Golestani R, Tio R, Zeebregts CJ, et al. Abdominal aortic calcification detected by dual X-ray absorptiometry: A strong predictor for cardiovascular events. *Ann Med.* 2010;42(7):539-545.
39. Lewis JR, Eggermont CJ, Schousboe JT, et al. Association between abdominal aortic calcification, bone mineral density and fracture in older women. *Journal of bone and mineral research : the official journal of the American Society for Bone and Mineral Research.* 2019.
40. Centers for Medicare and Medicaid Services. Physician Fee Schedule Look-Up Tool. 2019; <https://www.cms.gov/Medicare/Medicare-Fee-for-Service-Payment/PFSlookup/index.html>. Accessed July 23, 2019, 2019.
41. Szulc P, Blackwell T, Kiel DP, et al. Abdominal aortic calcification and risk of fracture among older women - The SOF study. *Bone.* 2015;81:16-23.
42. Szulc P, Blackwell T, Schousboe JT, et al. High hip fracture risk in men with severe aortic calcification: MrOS study. *J Bone Miner Res.* 2014;29(4):968-975.
43. Bagger YZ, Rasmussen HB, Alexandersen P, Werge T, Christiansen C, Tanko LB. Links between cardiovascular disease and osteoporosis in postmenopausal women: serum lipids or atherosclerosis per se? *Osteoporos Int.* 2006.
44. Flipon E, Liabeuf S, Fardellone P, et al. Is vascular calcification associated with bone mineral density and osteoporotic fractures in ambulatory, elderly women? *Osteoporos Int.* 2012;23(5):1533-1539.
45. Samelson EJ, Cupples LA, Broe KE, Hannan MT, O'Donnell CJ, Kiel DP. Vascular calcification in middle age and long-term risk of hip fracture: the Framingham Study. *J Bone Miner Res.* 2007;22(9):1449-1454.
46. Rodriguez AJ, Leow K, Szulc P, et al. Abdominal aortic calcification, bone mineral density and fractures: a systematic review and meta-analysis protocol. *BMJ Open.* 2019;9(4):e026232.
47. Schousboe JT, Richter SA, Beran MS. Potential Clinical Impact of Abdominal Aortic Calcification on Bone Density Lateral Spine Images. *J Clin Densitom.* 2016;19(4):436-443.

United States / Latin America
 250 Campus Drive
 Marlborough, MA 01752 USA
 Tel: +1.508.263.2900
 Sales: +1.781.999.7453
 Fax: +1.781.280.0668
www.hologic.com

MED-00410 Rev.001 (04/20) © 2020 Hologic, Inc. All rights reserved. Printed in USA. Specifications are subject to change without prior notice. Hologic, HD Instant Vertebral Fracture is a trademark and/or registered trademarks of Hologic and/or Hologic subsidiaries in the United States and/or other countries. All other trademarks and registered trademarks are the property of their respective owners. This information is intended for medical professionals in the U.S. and other markets and is not intended as a product solicitation or promotion where such activities are prohibited. Views and opinions expressed herein by third parties are theirs alone and do not necessarily reflect those of Hologic. Because Hologic materials are distributed through websites, eBroadcasts and tradeshow, it is not always possible to control where such materials appear. For specific information on what products are available for sale in a particular country, please contact your local Hologic representative or write to info@hologic.com